




Sample Name : M5861-Paroxetine HCl-02-DMSO
Acq. Operator : C22
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\2-POS-MON-1.M
Analysis Method : D:\METHODS\2-POS-MON-1.M
Sample Info : Easy-Access Method: '2-POS-MON-1'
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.5min, 95%B for 1.5 min
Flow Rate: 2ml/min
Column: Sunfire, 4.6*50mm, 3.5um
Oven Temperature: 50 C
MASS Range: 120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
0.52	0.02	12.30	10.67	0.57
1.30	0.02	7.49	6.59	0.35
1.45	0.02	2131.42	1995.98	98.48
2.52	0.02	13.14	8.81	0.61

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.45	0.02	215.35	215.07	100.00

Ret. Time: 1.47

<<< POSITIVE SPECTRA >>>

ES-API Positive

