

LCMS Report

Sample Information

LCMS

Acquired by : System Administrator
Sample Name : Genistein-02
Method File : B 5%-95%-P-new - 2.6min.lcm

Date Acquired : 2021/7/10 14:49:41


Injection Volume : 1

Vial : 104


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min
Pump A : 0.1% FA and 10% ACN in H2O
Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm


Peak#	Ret. Time	Height	Area	Area%
1	1.258	20467	41797	0.758
2	1.516	3534914	5475106	99.242
总计		3555381	5516902	100.000

PDA Ch2 214nm


Peak#	Ret. Time	Height	Area	Area%
1	1.255	16945	22742	0.653
2	1.516	2087802	3461572	99.347
总计		2104747	3484313	100.000

Mass Spectrum


Rt=0.169min 正




Rt=0.304min 正




Rt=0.462min 正




Rt=1.210min 负




Rt=1.480min 负




Rt=1.474min 正




Rt=2.110min 负



Rt=2.245min 负



Rt=2.402min 负



Rt=2.419min 正

