

LCMS Report

Sample Information

LCMS


Acquired by : System Administrator
Sample Name : Ibuprofen-02
Method File : B 5%-95%-P-new - 2.6min.lcm

Date Acquired : 2021/7/10 14:43:25
Injection Volume : 8
Vial : 102


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min
Pump A : 0.1% FA and 10% ACN in H2O
Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm

Peak#	Ret. Time	Height	Area	Area%
1	1.869	229469	397471	100.000
总计		229469	397471	100.000


PDA Ch2 214nm

Peak#	Ret. Time	Height	Area	Area%
1	1.872	2927470	7526947	100.000


Peak#	Ret. Time	Height	Area	Area%
总计		2927470	7526947	100.000

Mass Spectrum


Rt=1.497min 正




Rt=1.525min 负




Rt=1.587min 正




Rt=1.817min 负




Rt=1.834min 正




Rt=2.082min 正




Rt=2.110min 负



Rt=2.217min 正



Rt=2.245min 负



Rt=2.425min 负

