

Shimadzu LCMS-029 Report

Sample Information

LCMS


Acquired by : System Administrator
Sample Name : M3969-Polydatin-02
Method File : B 5%-95%-P-new - 2.6min.lcm

Date Acquired : 2021/11/3 13:18:56
Injection Volume : 2
Vial : 78


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min
Pump A : 0.1% FA and 10% ACN in H2O
Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm


Peak#	Ret. Time	Height	Area	Area%
1	1.208	1305991	2127182	98.766
2	1.385	16106	26572	1.234
总计		1322097	2153754	100.000

PDA Ch2 214nm


Peak#	Ret. Time	Height	Area	Area%
1	1.212	2721659	5896724	100.000
总计		2721659	5896724	100.000

Mass Spectrum


Rt=1.165min 负




Rt=1.182min 正




Rt=1.345min 负




Rt=1.362min 正




Rt=1.727min 负




Rt=1.744min 正




Rt=2.059min 正



Rt=2.087min 负



Rt=2.194min 正



Rt=2.267min 负

