



Agilent Technologies

Sample name: M5557-Deferasirox-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1E, Sig=240.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.715	MM m	0.09	1.15	0.39	0.01
7.984	MM m	0.14	3.09	0.75	0.03
9.736	MM m	0.74	11899.10	2356.56	99.96