

LCMS Report

Sample Information

LCMS

Acquired by : System Administrator

Sample Name : M9828-2-D08-01

Method File : B 5%-95%-P-new - 2.6min.lcm

Injection Volume : 15

Vial : 93


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min

Pump A : 0.1% FA and 10% ACN in H2O


Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm


Peak#	Ret. Time	Height	Area	Area%
1	1.382	1292738	1891966	99.426
2	1.522	12490	10929	0.574
总计		1305228	1902895	100.000

PDA Ch2 214nm


Peak#	Ret. Time	Height	Area	Area%
1	1.382	2632344	3892221	100.000
总计		2632344	3892221	100.000

Mass Spectrum


Rt=1.345min 负




Rt=1.339min 正




Rt=1.474min 正




Rt=1.632min 正




Rt=1.789min 正




Rt=1.975min 负




Rt=2.132min 负



Rt=2.239min 正



Rt=2.245min 负



Rt=2.425min 负

