

File D:\DATA\PRO3013-LCMS

Sample Info : Easy-Access Method: '1-POS-3MIN'


Method Info : Mobile Phase : A:Water(10mM NH₄HC0₃), B: Acetonitrile
Gradient : 5%B increase to 95%B within1.4min,
95%B for 1.6min,back to 5%B within 0.01min.

Flow Rate : 1.8mL/min


Column : Xbridge, 3.5um, 50*4.6mm

Oven Temperature : 50 C MS:103-1000


ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
0.51	0.01	10.63	18.56	0.77
0.53	0.02	15.05	10.14	1.09
1.28	0.02	1335.69	1059.15	98.14

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
0.51	0.01	2.01	5.07	0.17
1.28	0.02	1206.53	960.25	99.83

Ret. Time: 1.29


<<<< POSITIVE SPECTRA >>>>




Ret. Time: 1.95



Ret. Time: 2.06



Ret. Time: 2.19



Ret. Time: 2.23

