



11.31

8.883
8.868
8.406
7.917
7.911
7.895
7.888
7.574
7.555
7.525
7.367
7.347
6.894
6.879

4.008
3.874

2.815
2.502

0.969



Current Data Parameters
NAME AMG-458-01
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date 20160529
Time 10.38
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8278.146 Hz
FIDRES 0.126314 Hz
AQ 3.9584243 sec
RG 322.5
DW 60.400 usec
DE 8.00 usec
TE 298.2 K
D1 1.00000000 sec
TD0 1

==== CHANNEL f1 =====
NUC1 1H
P1 10.70 usec
PL1 -3.00 dB
SF01 400.1324710 MHz

F2 - Processing parameters
SI 32768
SF 400.1300032 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

