




Current Data Parameters
NAME AS703026-01
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date 20160531
Time 11.35
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8278.146 Hz
FIDRES 0.126314 Hz
AQ 3.9584243 sec
RG 456.1
DW 60.400 usec
DE 8.00 usec
TE 297.5 K
D1 1.00000000 sec
TD0 1

==== CHANNEL f1 =====
NUC1 1H
P1 10.70 usec
PL1 -3.00 dB
SF01 400.1324710 MHz

F2 - Processing parameters
SI 32768
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



9.375
8.508
8.210
8.198
7.713
7.709
7.683
7.670
7.491
7.303
7.281

3.966
3.683
3.669
3.656
3.651
3.638
3.625
3.444
3.430
3.418
3.410
3.398
3.384
3.369
3.365
3.356
3.350
3.337
3.205
3.190
3.172
3.157
3.139
2.510

0.007

1.000
1.029
1.067
1.045
2.099
1.062
1.057

1.339
3.634
1.195

