



File D:\DATA\
Sample Name : AZ 960-01
Acq. Operator : Inj : 0
Spec. Reported : MS Integration Inj Volume : 2 ul
Acq. Method : D:HPLC
Analysis Method : D:HPLC
Sample Info : Easy-Access Method:
Method Info : Mobile Phase: A: Water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%B increase to 95%B within 8min,
95%B for 7min, back to 5%B within 0.01min.
Flow Rate: 1 mL/min
Column: SunFire C18, 4.6*150mm, 3.5um
Column Temperature: 40 C

DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
8.49	0.06	2974.06	784.25	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
8.49	0.06	2146.84	567.94	100.00