



Current Data Parameters

NAME AZD2281-10
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters

Date_ 20180919
Time_ 11.11
INSTRUM spect
PROBHD 5 mm PADUL 13C
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.122266 Hz
AQ 4.0894465 sec
RG 199.46
DW 62.400 usec
DE 6.50 usec
TE 296.8 K
D1 1.00000000 sec
TD0 1

==== CHANNEL f1 =====

SFO1 400.1324710 MHz
NUC1 1H
P1 10.00 usec
PLW1 20.00000000 W

F2 - Processing parameters

SI 65536
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

1.06

1.00

1.02

2.06

2.05

1.02

2.07

5.16

1.12

2.16

1.07

0.35

4.08