

File D:\DATA\

Sample Name : AZD6244-06

Spec. Reported : MS Integration

Inj Volume : 2 ul

Acq. Method : D: HPLC

Analysis Method : D: HPLC

Sample Info : Easy-Access Method:

Method Info : Mobile Phase: A: water (0.01%TFA) B: ACN (0.01%TFA)

Gradient

Time 0 8 15


B % 5 95 95

Flow Rate: 1.0mL/min


Oven temperature: 40C

Column: SUNFIRE C18, 150*4.6mm,3.5um

DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
6.56	0.06	21.66	5.20	0.88
6.82	0.04	2429.62	938.41	98.34
7.79	0.04	19.40	6.91	0.79

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
6.82	0.04	1244.46	484.36	100.00