



LC/MS Report

Injection Date : 8 Apr 21
 Sample Name : AZD7762-02
 Method Info : Mobile Phase : A:water(0.01%TFA) B:ACN(0.01%TFA)
 Gradient:5%B increase to 95%B within 1.5min,95%B for 1.5min
 Flow Rate :2.0ml/min
 Column :SunFire C18 3.5um, 4.6mm*50mm
 Column Temperature:50C
 Detection:UV(214nm,4nm) and MS(ESI,POS Mode,110-1300amu)


*DAD1 A, Sig=2144 Ref=off



*DAD1 B, Sig=2544 Ref=off



MSD1 TIC, MS File ES-AP I, Pos, Scan, Frag: 100



Integration Results for DAD1 A, Sig=2144 Ref=off

RetTim	Width	Area	Height	Area%
1.16	0.02	589.09	516.33	100.00

Integration Results for DAD1 B, Sig=2544 Ref=off

RetTim	Width	Area	Height	Area%
1.16	0.02	827.91	725.16	100.00

Ret. Time: 1.17

<<<< POSITIVE SPECTRA >>>>

