

MASS SPECTROMETRY REPORT

X-Y-0626 4 (0.074) Cm (4-2)

Scan ES+
4.29e6



Sample Description

Analyzed date: 2019-06-27
Analyst: PY
Sample: Angiotensin II human-02
M.W.: 1046.2
Lot. No.: P181120-1

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Probe Bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN