


7.266  
5.620  
5.609  
5.599  
5.374  
5.359  
5.348  
5.340  
5.337  
5.326  
5.313  
5.298  
4.558  
4.541  
4.532  
4.526  
4.511  
4.502  
4.144  
4.129  
4.113  
4.099  
3.811  
3.790  
3.776  
3.754  
3.500  
3.473  
2.369  
2.361  
2.350  
2.342  
2.330  
2.311  
2.303  
2.016  
2.002  
1.989  
1.849  
1.628  
1.611  
1.593  
1.575  
1.300  
1.267  
0.897  
0.880  
0.863  
-0.000  
-0.008

Sample Name: DOTAP chloride-03  
Chemical Formula:  $C_{42}H_{80}ClNO_4$   
Molecular Weight: 698.56



CAS: 132172-61-3  
CDCl<sub>3</sub>, 400 MHz

