


8.141
8.138
8.121
8.119
7.876
7.873
7.855
7.838
7.835
7.725
7.705
7.567
7.548
7.530
7.472
7.452
7.444
7.441
7.433
7.425
7.303
7.280
7.272
7.183
7.163
7.075
7.057
7.038
6.912
6.905
6.882
6.873
4.751
4.080
4.073
4.062
4.055
4.038
4.022
4.007
3.989
3.983
3.972
3.965
3.741
3.724
3.317
3.188
2.506
2.501
2.497
2.465
2.445
2.393
2.378
2.365
1.987
1.976
1.249
1.232
1.123
1.106
1.089
-0.000

Sample Name:Erastin-03

Chemical Formula: C₃₀H₃₁ClN₄O₄

Molecular Weight: 547.05



CAS: 571203-78-6

DMSO-d₆, 400 MHz

1.00 1.02 1.00 2.08 1.04 2.04 2.03 1.04

2.01

2.07

1.01

4.09

2.14

2.13

3.11

3.13

10

8

6

4

2

0

PPM