



File D:\DATA\
 Sample Name : KPT-330-04
 Acq. Method : D: HPLC
 Analysis Method : D: HPLC
 Sample Info : Easy-Access Method:
 Method Info : Mobile Phase: A: Water(0.01%TFA) B: ACN(0.01%TFA)
 Gradient: 5%B increase to 95%B within 8min,
 95%B for 7min, back to 5%B within 0.01min.
 Flow Rate: 1 mL/min
 Column: SunFire C18, 4.6*150mm, 3.5um
 Column Temperature: 40 C



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
5.93	0.04	154.91	56.66	1.77
6.32	0.05	8581.56	2771.92	98.23

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
5.93	0.04	177.75	65.38	1.63
6.32	0.05	10752.24	3261.79	98.37