



Current Data Parameters
NAME KU-0063794-01
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date 20160529
Time 9.48
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8278.146 Hz
FIDRES 0.126314 Hz
AQ 3.9584243 sec
RG 362
DW 60.400 usec
DE 8.00 usec
TE 298.3 K
D1 1.00000000 sec
TDO 1

==== CHANNEL f1 =====
NUC1 1H
P1 10.70 usec
PL1 -3.00 dB
SF01 400.1324710 MHz

F2 - Processing parameters
SI 32768
SF 400.1300030 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



8.307
8.302
8.203
8.181
8.073
8.068
8.052
8.046
7.604
7.583
7.099
7.077

5.209

4.684
4.653
4.574

3.865
3.784
3.772
3.686
3.675
3.570
3.323

2.594
2.567
2.561
2.502

1.193
1.178

1.069
1.084
1.012
1.021
1.000

0.937
2.016
2.075

3.249
4.272
4.323
2.259

2.012

6.976

8 7 6 5 4 3 2 1 0 ppm