




Agilent Technologies

Sample name: M10048-PF05231023-02
Instrument: 2
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1D,Sig=235,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
5.508	MM m	0.04	4.53	1.84	0.08	
6.036	MB m	0.05	5464.19	1766.21	98.04	
6.316	MM m	0.04	1.27	0.56	0.02	
6.454	MM m	0.03	1.25	0.67	0.02	
6.624	MM m	0.05	102.34	35.13	1.84	
	Sum		5573.59			