



Sample name: M10495-AMI-1 disodium salt-01

Injection date:

Instrument: A010

Location:

Inj. volume: 2


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
4.275		0.3	0.2	0.002
4.349		3.2	1.0	0.022
4.473		11.3	4.4	0.077
4.565		14558.3	2239.5	99.600
5.569		27.7	6.2	0.190
5.759		8.2	2.6	0.056
6.112		1.6	0.5	0.011
6.404		4.2	1.0	0.029
6.622		1.8	0.6	0.012