




Sample name: M11540-Propranolol-01  
Instrument: 1  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
6.046	MM m	0.08	2.01	0.68	0.01
6.102	MM m	0.08	2.83	1.03	0.02
6.396	MM m	0.09	1.62	0.62	0.01
6.581	BB	0.42	14782.90	2884.47	99.56
7.123	MM m	0.08	1.18	0.42	0.01
8.876	MM m	0.23	44.23	8.45	0.30
8.920	MB m	0.09	13.84	4.92	0.09