



Sample name: M11545-Dichlorophen-01

Injection date:

Instrument: A010

Location:

Inj. volume: 2.5


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
7.626	MM m	0.05	72.33	30.66	0.71
9.831	BB	0.56	10030.65	2047.28	98.14
12.333	MM m	0.06	117.98	40.54	1.15