



Sample name: M13920-AZ084-01

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A,Sig=220.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.335		7.6	2.6	0.111
5.448		3.8	1.3	0.055
5.693		6728.5	2006.2	98.131
5.822		44.1	14.3	0.644
5.915		3.5	1.3	0.051
6.121		4.0	1.4	0.059
6.321		59.8	19.5	0.872
7.124		5.3	1.6	0.078