




Agilent Technologies

Sample name: M14076-Iptacopan-01
Instrument: 1
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
5.087	MB m	0.05	4.34	1.40	0.07	
5.834	MM m	0.07	3.35	0.86	0.05	
6.026	MM m	0.04	1.05	0.38	0.02	
6.199	MM m	0.05	0.80	0.28	0.01	
6.877	BB	0.30	6321.04	1977.93	99.50	
7.283	BB	0.13	6.63	2.01	0.10	
7.809	MB m	0.05	8.74	2.64	0.14	
8.051	MM m	0.05	6.64	2.22	0.10	