




Sample name: M14237-BAY-8002-01
Instrument: A010
Inj. volume: 1.5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
7.998		0.7	0.2	0.005
9.255		195.8	49.3	1.429
9.715		13438.8	2532.3	98.106
10.219		15.7	3.9	0.115
10.351		10.9	2.7	0.079
10.469		32.4	8.1	0.237
10.838		3.9	0.8	0.029