




Sample Name : M14920-Lenvatinb Mesylate-01
Acq. Operator : C22
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3MIN.M'
Method Info : Column:Xbridge C18(2) (4.6x 50 mm, 3.5um)
 Mobile phase: H2O(10 mmol NH4HCO3) (A) / ACN(B)
 Elution program: Gradient from 5 to 95% of B in1.5min at 1.8ml/min
 Temperature: 50 °C
 Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,103 to 800 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
0.81	0.02	29.96	18.06	1.61
1.18	0.02	9.33	7.88	0.50
1.36	0.02	1819.35	1163.53	97.89

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
0.81	0.02	23.73	14.29	1.25
1.18	0.02	8.69	7.07	0.46
1.36	0.02	1861.18	1175.46	98.29

Ret. Time: 1.37

<<<< POSITIVE SPECTRA >>>>

*ES-API Positive

