


MASS SPECTROMETRY REPORT



Sample Description

Analyst: SHAO
Sample: M15097- α -CGRP (mouse, rat) (TFA)-01
M.W.: 3806.3
Lot. No.:

Instrument: Waters 2695
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

A: 3806.95 \pm 1.01
Probe Bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.5ml/min
B. Conc.: 50%H₂O/50%ACN