




Sample name: M1787-INCB18424-05
Instrument: A010
Inj. volume: 1.5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1D, Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
5.752	MM m	0.08	3.59	1.40	0.06
5.812	MB m	0.09	8.03	2.70	0.13
6.360	MM m	0.29	5942.94	1751.21	99.46
7.202	MB m	0.11	3.72	1.29	0.06
7.817	MM m	0.14	16.71	5.29	0.28