

Shimadzu LCMS-029 Report

Sample Information

LCMS


Acquired by : System Administrator
 Sample Name : M1820-Tofacitinib-04
 Method File : B 5%-95%-P-new - 2.6min.lcm

Date Acquired : 2021/8/16 10:33:59
 injection Volume : 2
 Vial : 5


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min
 Pump A : 0.1% FA and 10% ACN in H2O
 Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm

Peak#	Ret. Time	Height	Area	Area%
1	1.093	173690	760186	100.000
总计		173690	760186	100.000


PDA Ch2 214nm

Peak#	Ret. Time	Height	Area	Area%
1	1.092	699270	3107465	100.000


Peak#	Ret. Time	Height	Area	Area%
总计		699270	3107465	100.000

Mass Spectrum


Rt=0.580min 负




Rt=1.075min 负




Rt=1.069min 正




Rt=1.497min 正




Rt=1.722min 正




Rt=2.087min 负




Rt=2.127min 正



Rt=2.245min 负



Rt=2.402min 负



Rt=2.419min 正

