




Agilent Technologies

Sample name: M1830-Firocoxib-01  
Instrument: 3  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
7.169		56.6	12.5	0.586
7.544	3.8	1.1	0.4	0.011
8.504	10.1	9595.4	2389.2	99.403