




Sample name: M1833-Imidaclopride-01  
Instrument: 2  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1F,Sig=260,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.544	BV	0.25	7.08	1.42	0.05
5.801	BM m	0.06	0.46	0.22	0.00
5.926	MM m	0.39	12795.94	3145.75	98.66
6.252	MM m	0.14	6.73	1.40	0.05
6.400	MM m	0.15	3.59	0.81	0.03
6.825	BB	0.27	145.71	39.73	1.12
7.125	BB	0.13	3.85	1.18	0.03
7.281	MM m	0.17	6.80	1.11	0.05