




Sample Name : M1985-Odanacatib-01-DMSO
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.5min,95%B for 1.5 min
Flow Rate:2ml/min
Column:Sunfire ,4.6*50mm,3.5um
Oven Temperature: 50 C
MASS Range:120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.50	0.03	17.18	9.55	0.93
1.63	0.02	10.84	7.10	0.59
2.08	0.01	1824.16	2016.36	98.49

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
2.08	0.01	2029.63	2315.93	100.00

Ret. Time: 0.38

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 2.10

