




Sample name: M20588-MD2-IN-1-01
Instrument: 1
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1B, Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
6.348	14.23	3.99	0.21
7.160	6.57	2.16	0.10
8.749	4.07	1.77	0.06
8.871	6863.08	1664.55	99.36
9.136	6.65	2.15	0.10
9.643	9.23	1.85	0.13
9.964	3.21	1.00	0.05