



Agilent Technologies

Sample name: M2218-Ribociclib (LEE011)-01

Injection date:

Instrument: 1

Location:

Inj. volume: HPLC


Type:

Acq. method:

Sample amount:

Processing method:

Manually modified:



Signal: DAD1C, Sig=254.0,4.0 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
5.487		0.24	2587.9	893.5	99.24
5.663		0.04	1.1	0.4	0.04
5.894		0.04	0.5	0.2	0.02
5.998		0.05	0.5	0.2	0.02
6.173		0.05	1.7	0.6	0.07
6.303		0.05	0.4	0.1	0.02
6.424		0.15	14.4	4.9	0.55
6.596		0.05	0.8	0.2	0.03
7.126		0.04	0.4	0.1	0.02
		Sum	2607.68		