



Sample name: M2225-Thalidomide-01

Instrument: A012

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
3.838	MM m	0.15	4.22	0.89	0.03
3.977	BB	0.13	5.58	1.54	0.04
4.781	MM m	0.13	5.90	1.70	0.04
5.428	MM m	0.10	1.98	0.67	0.01
5.979	BV	0.37	13586.51	2861.68	99.81
7.474	MB m	0.10	3.03	1.01	0.02
8.028	MM m	0.11	5.13	1.58	0.04