



Sample name: M22361-FTBMT-01

Instrument: A010

Inj. volume: 1.5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B,Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
6.594	57.33	16.59	0.68
7.068	2.75	0.88	0.03
7.177	9.28	3.30	0.11
7.508	2.42	0.62	0.03
7.866	8399.66	2315.14	99.09
8.242	5.48	1.05	0.06