


Sample Name : M2310-BX-795-03
Acq. Operator : 9
Spec. Reported : Integration

Location : P1-D-04
Inj : 0
Inj Volume : 5 ul

DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
2.84	0.05	19.70	6.57	0.26
3.51	0.05	37.94	12.83	0.51
3.75	0.07	7395.26	1585.57	98.65
4.21	0.06	14.33	3.26	0.19
4.50	0.06	29.28	6.63	0.39