




Sample name: M2352-10058-F4-03  
Instrument: 1.5  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1E,Sig=280,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
6.531	MM m	0.05	2.05	0.60	0.05	
8.687	MM m	0.07	3.34	0.79	0.08	
9.621	MM m	0.06	1.81	0.45	0.04	
10.603	MB m	0.07	10.96	2.33	0.27	
11.331	VM m	0.10	4055.45	585.68	99.55	
	Sum		4073.60			