




Sample name: M2748-GW9662-05  
Instrument: A011  
Inj. volume: 1  
Acq. method: HPLC  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1F,Sig=260,4 Ref=off

| RT [min] | Type | Width [min] | Area    | Height  | Area% |
|----------|------|-------------|---------|---------|-------|
| 5.747    | MM m | 0.16        | 0.54    | 0.12    | 0.01  |
| 6.314    | MM m | 0.15        | 2.32    | 0.61    | 0.02  |
| 6.754    | MM m | 0.15        | 0.50    | 0.14    | 0.01  |
| 7.029    | MM m | 0.10        | 0.50    | 0.17    | 0.01  |
| 8.891    | BB   | 0.45        | 9859.66 | 2434.43 | 99.96 |