




Sample name: M2996-Salubrinal-02
Instrument: 2
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1B,Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
4.533	10.34	3.79	0.23
5.535	5.79	1.97	0.13
7.108	5.98	1.30	0.13
7.288	1.87	0.75	0.04
7.408	2.90	1.02	0.06
8.755	4562.27	934.37	99.41
Sum	4589.16		