




Sample name: M3044-T0070907-03
Instrument: A012
Inj. volume: 2
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1F,Sig=260,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
3.952	MM m	0.10	0.26	0.09	0.00
4.224	MM m	0.13	1.16	0.32	0.01
5.117	MM m	0.10	0.45	0.16	0.00
5.667	MM m	0.47	13464.45	3187.44	99.97
6.985	MM m	0.13	1.55	0.45	0.01