



Sample name: M30897-Reproxalap-01

Injection date:

Instrument: A010

Location:

Inj. volume: 2.5


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1C,Sig=254.0,4.0 Ref=off

RT [min]	Peak Resolution	Width [min]	Area	Height	Area%
6.344		1.14	15738.2	3229.5	99.69
7.469		0.22	1.6	0.3	0.01
8.312		0.23	2.5	0.5	0.02
8.699		0.24	4.8	0.9	0.03
10.811		0.51	37.0	6.3	0.23
11.912		0.30	2.5	0.4	0.02