




Sample name: M3188-Istradefylline-01
Instrument: 1
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.865		1.8	0.6	0.045
6.040		0.9	0.3	0.023
6.168		4.9	1.7	0.126
8.918		3893.9	713.2	99.556
9.270		9.8	2.5	0.250
	Sum	3911.25		