




Sample name: M3189-TCN-201-01
Instrument: 2
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
7.119	MM m	0.06	2.56	0.63	0.03	
7.768	MM m	0.05	7628.28	2296.40	99.91	
8.160	MM m	0.06	2.67	0.67	0.04	
8.527	MM m	0.05	1.62	0.47	0.02	
		Sum	7635.13			