




Sample name: M3208-Geniposide-01
Instrument: 2
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1E, Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
4.110	BB	0.23	19.40	3.59	0.71	
4.381	MM m	0.06	1.44	0.39	0.05	
4.564	MM m	0.05	2727.59	852.48	99.13	
4.702	MM m	0.04	2.25	0.94	0.08	
4.762	MM m	0.05	0.65	0.23	0.02	
5.107	MM m	0.04	0.26	0.11	0.01	
		Sum	2751.60			