



Sample name: M3325-Captopril-02

Instrument: A010

Inj. volume: 5

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A, Sig=220,4 Ref=off

RT [min]	Width [min]	Area	Height	Area%
5.357	0.48	13595.8	1711.4	99.480
5.539	0.14	9.3	2.0	0.068
5.739	0.11	0.9	0.3	0.007
6.049	0.24	60.9	14.9	0.446