




Sample Name : M3649-GNE-7915-02
Acq. Operator : 11
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\2-POS-MON-1.M
Analysis Method : D:\METHODS\2-POS-MON-1.M
Sample Info : Easy-Access Method: '2-POS-MON-1'
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.5min, 95%B for 1.5 min
Flow Rate: 2ml/min
Column: Sunfire, 4.6*50mm, 3.5um
Oven Temperature: 50 C
MASS Range: 120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.96	0.01	1146.53	1222.80	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.96	0.01	1030.28	1096.77	100.00

Ret. Time: 1.98

<<<< POSITIVE SPECTRA >>>>

ES-API Positive

