



Sample name: M4572-Limonin-02

Instrument: A010

Inj. volume: 3

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1B, Sig=210.0,4.0 Ref=off

RT [min]	Area	Height	Area%
6.913	7.80	2.52	0.15
7.817	2.71	0.81	0.05
7.968	4.04	1.36	0.08
8.285	5098.01	1193.97	98.98
8.370	29.24	12.29	0.57
9.246	8.68	2.40	0.17