


Sample Name : M5212-GSK2879552-03- Methanol
Method Info : Mobile Phase: A:Water(10mM NH₄HC₃) B:ACN
Gradient: 5%B increase to 95%B within 8min,95%B for 7min
Flow Rate: 1.0ml/min
Column:Xbridge C18,3.5um,4.6*150mm
Column Temperature:45C
Detection:UV(214,4nm) and MS(ESI,Pos mode,110 to 1500 amu)



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
5.55	0.05	47.81	13.74	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
5.55	0.04	4322.44	1813.62	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
5.55	0.04	574.96	246.66	100.00