




Sample name: M54520-TC11-01
Instrument: A010
Inj. volume: 2.5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1C, Sig=254.0, 4.0 Ref=off

RT [min]	Peak	Width [min]	Area	Height	Area%
10.249		0.22	3.1	0.5	0.02
10.407		0.16	1.2	0.2	0.01
10.810		0.73	16813.5	3054.4	99.48
11.396		0.32	16.4	1.5	0.10
12.052		0.47	66.5	12.1	0.39